

## Helmholtz Coils for MRP

The length of the saturated absorption cell is  $l_c = 8$  cm and has a radius of  $r_c = 1.25$  cm. In order to create a magnetic field around this, we will use Helmholtz coils. The magnetic field of Helmholtz coils is given by

$$B = \left(\frac{4}{5}\right)^{3/2} \frac{\mu_0 n I}{R}. \quad (1)$$

We want to have the following conditions satisfied:

- $I = 1$  A
- $0 \text{ g} \leq B \leq 10 \text{ g}$
- $R = 27/2 \text{ cm} = 0.135 \text{ m}$

Solving for the number of turns in each Helmholtz coil  $n$  to create the maximum magnetic field of 10 gauss, we find

$$\begin{aligned} n &= \left(\frac{5}{4}\right)^{3/2} \frac{RB}{\mu_0 I} \\ &= \left(\frac{5}{4}\right)^{3/2} \frac{(0.135 \text{ m})(0.001 \text{ T})}{(4\pi \times 10^{-7} \text{ Tm/A})(1 \text{ A})} \\ &\approx 150 \text{ turns} \end{aligned} \quad (2)$$

With of wire gauge of about 20 AWG (is this about standard?), the wire will have a diameter of about  $w_d = 1$  mm and will be able to hold up to 10 A. The width of the coil will then need to be

$$\begin{aligned} L_{coil} &= \left(\frac{n}{x}\right) w_d \\ &= \left(\frac{150}{3}\right) 1 \text{ mm} \\ &= 5 \text{ cm} \end{aligned} \quad (3)$$

where  $n$  is the number of turns and  $x$  is the number of times the wire is layered (3 layers in this situation).

We will thus need length of wire

$$\begin{aligned} l_{wire} &= n * \text{circumference} * 2 \text{ coils} \\ &= 150 \text{ turns} * 2\pi (27/2 \text{ cm}) * 2 \\ &= 250 \text{ m} \end{aligned} \quad (4)$$